

EXPERIENCE THE HEALING NATURE OF LIGHT

med^{HOME}



NEW MEDX^{HOME} 1000mW LIGHT THERAPY TREATMENT SYSTEM

The MedX^{HOME} uses high power medical near infrared (870nm) and red visible (633nm) LEDs. The MedX^{HOME} unit has been developed to offer a home care choice to patients.

This unit is lightweight and user-friendly and with hands-free superluminous diode (SLD) technology to administer light therapy treatments in the comfort of the patient's home. Protocols are programmed by the health professional.

A maximum of five different treatments can be programmed into the unit at one time. The number of days for use is set by the clinician to ensure it and the patient return to the clinic. The unit set-up is quick and easy and makes preparing for treatments simple.

MedX^{HOME} as a supplement to the MedX Laser and Light Therapy treatment System is an additional revenue stream with the return exceeding the original investment within the first month of renting the units to patients. The unit comes in a practical carrying case with supporting educational documentation for both clinician and patient.


MedX products are FDA cleared and Health Canada Licensed under following numbers: 15556, 76086, 69680 & 60952. Our products are produced in an ISO 13485:2003 certified manufacturing facility in Canada.

MedX Light Therapy is a proven and highly effective treatment which increases local microcirculation, accelerates tissue repair, reduces pain, and decreases inflammation/edema.



1-888-363-3112 (US/Cdn only) • 1-905-670-4428 • info@medXhealth.com • www.medXhealth.com

EXPERIENCE THE HEALING NATURE OF LIGHT



MEDX^{HOME} LIGHT THERAPY TREATMENT SYSTEM

INCLUDES:

MedX^{HOME} unit:

A compact, light-weight unit with one SLD connection port.
 LCD screen with easy to read instructions and 4 operating buttons.

SLD cluster:

Superluminous diode (SLD) cluster provides a total of 1000mW of energy.
 SLD covers large surface area of 22.48cm².

With the included 18" flexstrap, the MedX cluster can be strapped to the body and used hands-free, or comfortably handheld over the area treated.

Standard components:

- Programming Key
- MedX^{HOME} Unit
- SLD Cluster
- AC Power Adapter
- 18" Flexstrap
- Carrying Case

Clinician's Kit

- 1 Programming & Operating Guide
- 1 Quick Reference Sheet
- 1 Patient Operating Instructions Sheet
- 1 Sample Rental Agreement Form
- 1 Warranty/ Registration Card
- 5 Home Instruction Sheets and Daily Home Record Sheets

MedX Kits: Product Description

| | |
|------|-------------------------------------|
| 0915 | MedX 450mW Portable Laser |
| 2031 | MedX ^{HOME} System |
| 3116 | MedX Console + 450mW Laser + 2 SLDs |

©MedX Health Corp. 2005

